



CAD/CAM SYSTEMS | INSTRUMENTS | HYGIENE SYSTEMS | **TREATMENT CENTERS** | IMAGING SYSTEMS

HUGO – THE SMART STOOL FOR YOUR DENTAL PRACTICE

Excellent ergonomics.

The Dental Company

sirona.

HUGO – INNOVATIVE DESIGN FEATURES

Intelligent design for enhanced comfort.

Sophisticated design delivers optimum freedom of movement and outstanding comfort. Intuitive, aesthetic, comfortable and hygienic. Simply sit down and start treatment.



HUGO Manual:
Manual height adjustment.
“Vital” colour combination



HUGO Manual Plus:
Manual height adjustment.
Foot ring. “Natural” colour
combination



HUGO Freehand:
Height adjustment via
foot-release. “Elegant” colour
combination



HUGO – GROUND-BREAKING TECHNOLOGIES

Ideas that move you ahead.



Backrest – dynamic support

The flexibly mounted HUGO backrest follows your movements and provides continuous support – while you work and while you relax. Your spine is automatically straightened, thus eliminating the risk of muscle strain.



Seat cushion – enhanced comfort

Thanks to the innovative OptiMotion system, the inclination of the front left-hand and right-hand half of the seat cushion is freely adjustable. This increases your freedom of movement and ensures excellent access to the patient. At the same time, the pressure on your thighs is significantly reduced. The new cool upholstery delivers enhanced comfort during prolonged treatment sessions.




Stool base – increased flexibility

The HUGO stool can be adjusted to suit the individual size of each member of your practice team. The foot-controlled height adjustment function (Version Freehand) is designed for hygienic hands-free operation.



HUGO – OPTIMISED ERGONOMICS

Adapted to your requirements from head to foot.



- 1 The smart, ultra-comfortable backrest supports your spinal column from top to bottom in each successive treatment position.
- 2 A spin-off from the space industry, the high-tech seat upholstery dissipates body heat. This means that you stay cool even during prolonged treatment sessions.
- 3 The unique OptiMotion system supports your natural movements, eliminates pressure points, and increases your room for manoeuvre.
- 4 HUGO dispenses with knobs and handles and adapts itself automatically to each sitting position. This allows you to maintain

a healthy posture at all times and minimise physical strain. In addition, the stool design supports effective infection control procedures.

- 5 At no extra charge you can choose between two different gas springs. This enables you to select the ergonomically correct stool height, depending on how tall you are.


Height of user* (cm)	Seat height
150–180	Short (405–540 mm)
170–205	Long (505–645 mm)


** These figures can vary from user to user!*

- 6 The protective pads on the stool base can be replaced in the event of heavy wear and tear.
- 7 The foot ring (Version Manual Plus) caters for users of different sizes. You and your assistant can both maintain an ergonomically correct posture.
- 8 The new foot-release (Version Freehand) permits the hands-free and hygienic adjustment of the seat height.

The Sirona Ergonomics Program.

The high-tech HUGO stool is a key element in the Sirona Ergonomics Program. It enables you to maintain a healthy posture throughout the treatment process. Dynamic, comfortable and versatile, HUGO is the perfect complement to your Sirona treatment centre.





CAD/CAM SYSTEMS | INSTRUMENTS | HYGIENE SYSTEMS | TREATMENT CENTERS | IMAGING SYSTEMS

SIRONA – UNIQUE WORLDWIDE SYSTEMS EXPERTISE IN DENTAL EQUIPMENT PRODUCTS

Sirona develops and manufactures a comprehensive range of dental equipment, including CAD/CAM Systems for dental practices (CEREC) and laboratories (inLab), Instruments and Hygiene Systems, Treatment Centers and Imaging Systems. Sirona manufactures high technology products that guarantee ease of use and a high return on investment – for the good of your practice and for the benefit of your patients. In this way, you can approach every challenge that you face, confident in the knowledge that: **It will be a great day. With Sirona.**

Sirona Dental Systems · E-mail: contact@sirona.com · www.sirona.com

The Dental Company

sirona.