

INNOVATIVE LED TECHNOLOGY

The new LEDview – Viewing with Precision.





INNOVATIVE LED TECHNOLOGY

The new LEDview – Viewing with Precision.

- Bright and consistent illumination for the **highest levels of precision**
- Three-joint support arm's for **flexible positioning**
- **High-quality design** and workmanship



The new LEDview – additional benefits and technical data.

Benefits

- Daylight illumination with a light intensity of up to 30,000 lux
- Anti-glare for the patient through clearly confined luminous spot
- Hygienic, fully enclosed casing
- Long-life LEDs eliminate the need for bulb replacement
- High efficiency and low power consumption
- No cooling fan required
- Hands-free operation by sensor

Technical Data

- Brightness dimmable at several levels (analog SIROLUX F)
- Composite circuit with reduced light intensity
- Color temperature > 5,700 K
- LED life-time > 20,000 h
- Power consumption < 10 W
- Sensor switch at lamp crown ("On", "Off", "Composite")



CAD/CAM SYSTEMS | INSTRUMENTS | HYGIENE SYSTEMS | TREATMENT CENTERS | IMAGING SYSTEMS

SIRONA – UNIQUE WORLDWIDE SYSTEMS EXPERTISE IN DENTAL EQUIPMENT PRODUCTS

Sirona develops and manufactures a comprehensive range of dental equipment, including CAD/CAM Systems for dental practices (CEREC) and laboratories (inLab), Instruments and Hygiene Systems, Treatment Centers and Imaging Systems. Sirona manufactures high technology products that guarantee ease of use and a high return on investment – for the good of your practice and for the benefit of your patients. In this way, you can approach every challenge that you face, confident in the knowledge that: **Enjoy every day. With Sirona.**

Your Sirona dealer:

Sirona Dental Systems · E-mail: contact@sirona.com · www.sirona.com

The Dental Company

sirona.