


AVAILABLE KITS



- T4 LINE KIT BMT / BOT**
- Including T4 LINE B 40
 - Including T4 LINE BH 40
 - Including Air Motor
 - Including T4 Racer



Type	REF.
T4 LINE Kit BMT	64 18 318
T4 LINE Kit BOT	64 18 326



- T4 LINE KIT BM / BO**
- Including T4 LINE B 40
 - Including T4 LINE BH 40
 - Including Air Motor



Type	REF.
T4 LINE Kit BM	64 02 049
T4 LINE Kit BO	64 02 056



- T4 LINE KIT BI**
- Including T4 LINE B 40
 - Including T4 LINE BH 40

Type	REF.
T4 LINE Kit BI	64 02 064

All products in this brochure are supplied without a paper version of the instruction manual. This can be easily downloaded in pdf format from the internet at www.sirona.com/manuals.

ALWAYS AT THE FOREFRONT OF INNOVATION!

As global innovation leader for dental equipment, we continuously invest in research and thus in the future of modern dentistry. By networking digital technologies with integrated solutions and optimizing the treatment workflow, we create improved treatment results, more comfort and safety for the patient as well as time and cost savings in everyday work. The combination of constant innovative power and globally growing sales and service structures makes Sirona the global market leader trusted by thousands of practices and labs around the world. **Enjoy every day. With Sirona.**



CAD/CAM systems
From pioneer to new standard. For almost 30 years we have been developing digital dentistry and creating new possibilities for the future practice and lab.



Imaging systems
Best image quality with the lowest dose. More than 100 years of developing x-rays for the dental practice make us the number 1 innovation partner.



Treatment centers
The business card of modern practices. We are striving to create the ideal ergonomic and innovative center. Individually tailored to the well-being and demands of the patient and dentist.



Instruments
Advantages that speak for themselves. We make sure that we provide the right balance of proven quality, individual ergonomics and innovative technology for user-friendly work.



Hygiene systems
Competence that gives you safety. When it comes to hygiene in the practice, we do not take any shortcuts.

RTS-RIEGELTEAM.DE
Subject to technical changes and errors in the text, Order No. AG1.100-M45-8643-01-7600, Printed in Germany, Dispo-No. 04603, JP 13-0597 WS 02131.V5

T4 BASIC.

WE ONLY PUT THE BEST INTO YOUR HANDS.



SIRONA.COM

The Dental Company

sirona.

The Dental Company

sirona.

YOUR CHOICE, FOR ESSENTIAL WORKING.

The instruments of the T4 class are the entrance to the Sirona world.
A clever combination of proven technology and engineering excellence
for a cost-effective use.

Quality

MADE IN GERMANY

Essential and Ergonomic

DESIGN

Can be sterilized and heat desinfected

HYGIENE

T4

BASIC



STRAIGHT AND CONTRA- ANGLE HANDPIECES



T4 LINE – BASIC PRICE – RELIABLE QUALITY

The straight and contra-angle handpiece of the T4 class delivers with an essential design and a robust steel sleeve. Due to the air motor – available with Midwest 4-hole or Borden 2-hole – you will stay flexible, regardless of the treatment unit you are using.

- Essential design out of stainless steel
- External Spray
- Push button chuck
- For use with burs that have a 2.35 mm bur shaft diameter



TECHNICAL DATA OVERVIEW*

	T4 LINE B 40	T4 LINE BH 40	T4 Air Motor BM / B0
Speed (rpm)	max. 40,000 18.000 – 21.000 if connected to the T4 Air Motor BM / B0	max. 40,000	18,000–21,000 clockwise/counterclockwise
Spray	External	–	External
Head diameter (mm)	8.7	–	–
Head height (mm)	12.4	–	–
Light/illuminance (lux)	–	–	–
Weight (g)	72	77	90
Recommended burr size (ISO)	2.35 mm	2.35 mm	–
REF.	64 01 827	64 01 843	64 32 947 / 64 32 954
REF. Package of 25	64 02 106	64 02 114	64 32 962 / 64 32 970

* All technical specifications are approximate values.

TURBINES



T4 TURBINES – BASIC PRICE – RELIABLE QUALITY

The highspeed handpiece of the T4 class convinces with its combination of an essential design and exceptional high power level of up to 17 watts. Regardless which treatment unit you use, with Midwest 4-hole and Borden 2-hole you have always the right choice.

- Essential design in a chrome-plated version
- Ceramic ball bearings
- Single Spray
- Up to 17 watts
- Push button chuck



TECHNICAL DATA OVERVIEW*

	T4 Racer Midwest	T4 Racer Borden
Speed (rpm)	400.000	400.000
Spray	Single	Single
Head diameter (mm)	11.4 mm	11.4 mm
Head height (mm)	14.5 mm	14.5 mm
Light/illuminance (lux)	–	–
Weight (g)	43	43
Recommended burr size (ISO)	1.6 mm	1.6 mm
REF.	63 71 467	64 05 489
REF. Package of 25	63 86 234	64 12 816

* All technical specifications are approximate values.